












































# Abutment Torque Reference Sheet

**Have questions? Need parts? Contact Technic Dental Lab!!!**

**Front Desk: (708)479-1212    Email:Technicdental@sbcglobal.net**

<u>Company</u>	<u>Implant System</u>	<u>Torque Values</u> (Ncm)		<u>Driver Tip</u>
	BioHorizons	30 Ncm		.050/1.25mm Hex
	Gold Screw Titanium	35 Ncm 20 Ncm	 	Gold-Square Driver Titanium .048 Hex
	Prosthetic/Screw 3.3 3.8 - 6.0	15 Ncm 20 Ncm 30 Ncm	  	.050/1.25mm Hex .050/1.25mm Hex .050/1.25mm Hex
	Ankylos Astra 3.0 Astra 3.5/4.0 (Aqua) Astra 4.5/5.0 (Lilac) Astra E.V. Frident/Xive	15 Ncm 15 Ncm 20 Ncm 25 Ncm 25 Ncm 24 Ncm	     	1.0mm Hex .050/1.25mm Hex .050/1.25mm Hex .050/1.25mm Hex .050/1.25mm Hex .048/1.2mm Hex
	HiOssen (Osstem)	Mini - 20 Ncm Regular - 30 Ncm	 	.048/1.2mm Hex

<u>Company</u>	<u>Implant System</u>	<u>Torque Values</u> (Ncm)		<u>Driver Tip</u>
	Legacy Replant Swishplant Screw Plus	30 Ncm 35 Ncm 35 Ncm 30 Ncm	   	1.25mm Hex 1.25mm Hex or Uni-grip 1.25mm Hex or SCS Driver 1.25mm Hex or SCS Driver
	External Hex Genesis Prima Renova Restore	30 Ncm 30 Ncm 30 Ncm 30 Ncm 30 Ncm	 	Titanium .048 Hex/ Gold Screw-QuadDriver
	MegaGen	35 Ncm		.048/1.2mm Hex
	Gold Titanium	30 Ncm 30 Ncm	 	.050/1.25mm Hex .050/1.25mm Hex
	Abutment Screw Branemark, Sterioss Multi-Unit Screw, Prosthetic Screw Nobel 3.0	35 Ncm 35 Ncm 15 Ncm 15 Ncm	   	Uni-Grip Driver Uni-Grip Driver Uni-Grip Driver Uni-Grip Driver
	Abutment Prosthetic	35 Ncm 10 Ncm	 	SCS/Torx Driver SCS/Torx Driver
	Abutment Prosthetic	30 Ncm 20 Ncm	 	.050/1.25mm Hex .050/1.25mm Hex